

Technical Data

cardiofax 

ECG-3350

Electrocardiograph

(Order code: ECG-3350K, made in Japan)

Fighting Disease with Electronics



NIHON KOHDEN

TD. ECG-3350

This technical data may be revised or replaced by Nihon Kohden at any time without notice.
Some products may not be available in your country.

FEATURES

ECG-3350 is designed to ensure Quality diagnosis, Quality examination, and Quality accessibility.

Quality diagnosis

- Complies with the IEC60601-2-25 international standard for electrocardiographs.
- Nihon Kohden's advanced filtering technology eliminates noise with less waveform distortion.
- Nihon Kohden's ECG analysis program, ECAPS12C provides more than 200 findings with detailed analysis including Brugada type ECG interpretation.
- Has function of synthesized 18-lead ECG, synECi18 which is Nihon Kohden original technology and 6 synthesized leads, V3R, V4r, V5R, V7, V8 and V9 are derived from the standard 12-lead ECG waveforms. It can be helpful to find ischemia on right ventricular or posterior wall where have the difficulty to find it with standard 12-lead ECG.

Quality examination

- 8.0 inch color display realizes clear 12-lead ECG monitoring. Alert message and guide function help operators check cause of error. Electrode detachment and incorrect placement for chest and limb leads can be detected and indicated as the alert message.
- The most recent 3 minutes of all lead ECG waveforms can be stored internally, and a part which should be recorded and analyzed can be selected.

Quality accessibility

- Data can be output via wired LAN or built-in wireless LAN module.
- A USB memory or an SD memory card can be used as an external storage.
- Patient information can be received from external server. And report can be output in DICOM or PDF format which realize greater data integration with various data management system.

SPECIFICATIONS

◆ ECG Input

Input impedance: $\geq 50\text{M}\Omega$ at 0.67Hz
Polarization-proof voltage: $\pm 550\text{ mV}$
Input protection function: defibrillation-proof function
Common mode rejection ratio: $> 105\text{dB}$
Input circuit current: $< 0.05\text{ }\mu\text{A}$
Standard sensitivity: $10\text{ mm/mV} \pm 2\%$
Internal noise: $\leq 20\text{ }\mu\text{Vp-p}$
Interference between channels: $\leq -40\text{dB}$
Frequency characteristics: with 10 Hz as benchmark, 0.05 to 150 Hz ($+0.4\text{ dB} / -3.0\text{ dB}$)
Sample rate: 16,000 samples/s/ch

Heart rate count measurement:

Measurement and display range: 30 to 300bpm
Accuracy: $\pm 10\%$ (30 to 240bpm)

◆ Waveform Data Processing

Number of leads: 12 leads
Input channel: 1 channel
Sample rate: 500 samples/s, $1.25\text{ }\mu\text{V/LSB}$
Response to minimum signal: $\leq 20\text{ }\mu\text{Vp-p}$
EMG filter: 25 / 35 Hz (-3dB)
High-cut filter: 75, 100, 150 Hz (-3 dB)
AC line filter: 50 / 60 Hz ($\geq -20\text{dB}$)
Baseline drift suppression: Weak: $0.1\text{Hz}(-20\text{dB})$, Strong: $0.1\text{Hz}(-34\text{dB})$
Time constant: $\geq 3.2\text{ s}$
Sensitivity: 5, 10, 20 mm/mV
Pacemaker detection: Available
A/D Conversion: 24bit

◆ Display

Size: 8.0 inch
Type: color TFT LCD
Resolution: 800 x 480 dots
Displayed data: 12 lead ECG waveform, patient information, recording settings, operation mode, heart rate, QRS sync mark, error message, electrode detachment, noise

◆ Recorder

Recording mode: Thermal line
Printing density: 200 dpi
Scanning line density: 1 ms
Number of recording channels: 3, 3+Rhythm, 6, 12
Recording speed: 5, 10, 12.5, 25, 50 mm/s
Recording speed accuracy: $\leq \pm 5\%$
Recording data: ECG waveform, heart rate, lead name, version, date and time, paper speed, sensitivity, filter setting, patient information, measured information, and marks
Recording paper: 210 mm width, 30 m long Z fold / 216mm width,
Extended recording: Available
Grid Recording: Available

◆ External Input / Output

External input: 2 channels, 10 mm/0.5 V, input impedance $\geq 100\text{k}\Omega$
Signal output: 1 channel, $0.5\text{ V/mV} \pm 5\%$, output impedance $\leq 100\Omega$

◆ ECG Analysis

Program name: ECAPS 12C
Analysis patient age: neonate, child, adult
Finding items: approx.200
Judging items: 5

◆ Data management

Output file: DICOM, PDF(with SD card/USB memory)
External printer: Available

◆ Power Requirements

Voltage: AC (100 to 240 V) $\pm 10\%$
Frequency: (50 / 60 Hz) $\pm 2\%$
Power input: ≤ 150 VA
Battery pack: 12V
Battery operation time: more than 180 minutes
Battery charging time: within 3 hours

◆ Interface

USB 1.1 A type: 2
SD card slot: 1
LAN port: 1

◆ Environment

Storage and transport environment

Main Unit

Temperature: -20 to +65°C
 Relative humidity: 10 to 95%
 Atmospheric pressure: 700 to 1,060 hPa

Battery

Storage temperature range
 Storing for under 1 month: -20 to +55°C
 Storing for over 1 month: -20 to +35°C
 Storage humidity range
 Storing for under 2 months: 20 to 85%
 Storing for over 2 months: 45 to 85%

Recording paper

Storage temperature range: 0 to 45°C
 Storage humidity range: 40 to 80%

Operating Environment

Main unit

Temperature: 5 to 40°C
 Relative humidity: 25 to 95%, noncondensing
 Relative humidity: 40 to 80%, recording paper
 Atmospheric pressure: 700 to 1,060 hPa

Battery

Operating humidity range: 40 to 85%

◆ Performance standard

IEC60601-2-25:2011

◆ Electromagnetic Compatibility

IEC 60601-1-2: 2014

◆ Safety

Safety Standard IEC 60601-1 3.1st edition: 2012, IEC 60601-2-25: 2011
Type of protection against electric shock: AC power: CLASS I
 Battery pack power: Internally powered equipment
Degree of protection against electric shock: Defibrillator proof type CF applied part
Degree of protection against harmful ingress of water: IPX0 (non-protected)
Safety grade of use in the environment: Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide
Mode of operation: continuous operation

◆ Dimensions and Weight

Main unit

Dimensions: 255W x 420D x 140H mm (including handle)
 Weight: 4.7 kg

STANDARD ACCESSORIES

Recording paper, FQW210-2-140 1
 Thermal head cleaner pen, Y011 1
 Input check jig, 6114-092972 1

OPTIONAL ACCESSORIES

Accessory kit,

Description	Q'ty	Code No.	Supply Code No.	YD-121D	YD-122D	YD-123D	YD-124D
Patient cable							
BA-901D, AHA 3 mm diameter tip	1	706256	K019			✓	
BI-901D, IEC 3 mm diameter tip	1	687864	K079	✓	✓		
BI-902D, IEC/DIN 4 mm diameter tip	1	705248	K081				✓
Power cord							
Type H	1	186656	L935		✓	✓	
Type N	1	314839	L936	✓			✓
Chest electrode, 3 mm diameter tip, for adult, 3 pcs/set	2	709957	H041	✓	✓	✓	
Chest electrode, 4 mm diameter tip, for adult, 3 pcs/set	2	666949	H043				✓
Clip-on limb electrode, 3 mm diameter tip, 4 pcs/set, Fastclip	1	709939	H068	✓	✓	✓	
Clip-on limb electrode, 4 mm diameter tip, 4 pcs/set, Fastclip	1	709948	H068K				✓

Wireless LAN module, QI-330D

Cart, KD-107D (Height: 790mm) / **KD-107E** (Height: 890mm)

Patient cable hanger, KH-200D

WLAN converter holder, DI-021D

Magnetic card reader holder, DI-022D

Middle shelf, DI-023D

Drawer, DI-024D

Tray set, DI-025D

Roll paper holder, DI-026D

Barcode reader holder, DI-027D

Synthesized 18-lead program, QP-230EK

◆ Local Purchase Options

Purchase the following recommended models.

Recommended models

1D Bar code scanner

- LS2208 USBR (ZEBRA)

<https://www.zebra.com/ap/en/products/scanners/general-purpose-scanners/handheld/ls2208.html>

If you purchase the recommended model, it is necessary to purchase three items of below together.

Main body: LS2208-SR20111R

USB cable: CBA-U01-S07ZAR

Manual: JP-72-58808-01J1

- Voyager 1250g-USB (Honeywell)

<https://www.honeywellaidc.com/en/products/barcode-scanners/general-duty/voyager-1250g>

2D barcode reader

- Xenon 1900GSR-1USB (Honeywell)
<https://www.honeywellaidc.com/products/barcode-scanners/healthcare>

Magnetic card reader

- CRF-200U-5101-00 (NEURON)
http://neuron.co.jp/English/index_E.htm

USB flash drive

- DTSE9G2/32G (Kingston)
<https://www.kingston.com/en/usb-flash-drives/datatraveler-se9-g2-3-0-usb-flash-drive>

- SDCZ50-016G (SanDisk)
<https://shop.westerndigital.com/products/usb-flash-drives/sandisk-cruzer-blade-usb-2-0#SDCZ50-128G-A46>

Reference specifications*

Bar code scanner

USB1.1, bus powered, HID class type

Magnetic card reader

USB1.1, bus powered, HID class type

*If the recommended model is not available, you can try using a unit with these reference specifications although the proper function of any unit other than the recommended models cannot be guaranteed.

CONSUMABLES

Battery:

[X065] Rechargeable battery, NKB-301V

Memory card:

[Y154D] SD memory card, QM-001D

Recording paper:

[A052] FQW210-10-295, 5 pcs/set, Z-fold, 210mm x 100m (with KD-107D/E+DI-023D)

[A053] FQW210-30-295, Z-fold, 210mm x 300m (with KD-107D/E+DI-023D)

[A054] FQW210-3-140, 10 pcs/set, Z-fold, 210mm x 30m

[A134K] Recording paper, RQS210-6, 10 pcs/set, roll, 210 mm x 60m (with KD-107D/E + DI-026D)

Patient cable:

[K111] Patient cable, BJ-961D, IEC, 3mm dia

[K079] Patient cable, BJ-901D, IEC, 3mm dia

[K081] Patient cable, BJ-902D, IEC, 4mm dia

[K082] Patient cable, BJ-903D, IEC, clip

[K019] Patient cable, BA-901D, AHA, 3mm dia

[K020] Patient cable, BA-903D, AHA, clip

Clip-on limb electrode

[H068A] 4 pcs/set, adult, for 3mm dia tip, FASTCLIP

[H068B] 4 pcs/set, adult, for 4mm dia tip, FASTCLIP

[H068D] 4 pcs/set, child, for 3mm dia tip, FASTCLIP

[H068E] 4 pcs/set, child, for 4mm dia tip, FASTCLIP

Chest electrode

[H041A] 3 pcs/set, adult, for 3mm dia tip, gray, suction

[H043A] 3 pcs/set, adult, for 4mm dia tip, suction

[H042C] NC-123D, 3 pcs/set, child, for 3mm dia tip, suction Octopus

[H042D] NC-123D, 3 pcs/set, child, for 3mm dia tip, suction Bear

[H044A] NC-124D, 3 pcs/set, child, for 4mm dia tip, suction Octopus

[H044B] NC-124D, 3 pcs/set, child, for 4mm dia tip, suction Bear

[H042E] NC-123D, child, for 3mm dia tip, 6 colors

[H044C] NC-124D, child, for 4mm dia tip, 6 colors

ECG past

[F010] Z-101BC, 100g tube x 2/set, CARDIOCREAM

Disposable Electrodes Vitrode

[G207] L-150X, 30 x 5 pcs/set, adult, 35mm dia, X-ray

[G270A] NC-032Y, 100 x 5/set, 20 x 22.5mm, tab type

ORDERING INFORMATION

Order code of the main unit: ECG-3350K
Country of origin: Japan

When ordering, please specify the following information:

Electrocardiograph model

Requirements for consumables

Requirements for optional accessories

Specify the code in [xxxx]



NIHON KOHDEN CORPORATION

1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8000 Fax +81 (3) 5996-8100
www.nihonkohden.com